


Stabilization of Mouse Embryonic Stem Cells



<p>Cultivation</p>	<p>15% FBS/DMEM (1 mM of Sodium pyruvate, 100 µM of NEAA, 100 µM of β-ME, 1000 U/ml of LIF) was used as culture medium. Mouse Embryonic Fibroblasts (MEF) were used as „feeder cells“.</p>
<p>Freezing</p>	<p>Cells were frozen in 5 vials / (60mm dish corresponds to 3.0×10^6 cells/vial). 1 ml/vial of Bambanker™ freezing medium was added and the mixture was directly frozen in -80 °C. The following day, the vials were transferred to liquid nitrogen (slow freezing).</p>
<p>Thawing</p>	<p>Cells were incubated at 37 °C and transferred in cooled culture media. After collection, cells were seeded in 6 well plates and 6 cm dishes.</p>
<p>Results</p>	<p>Stabilization of Mouse Embryonic Stem Cells by using Bambanker™ was succesful. Cells were undifferentiated, even after freeze and thaw procedure. No modifications of cells could be observed. <i>Data were kindly provided by Dr. Ahn (Tokyo Institute of Technology Graduate School of Bioscience and Biotechnology Department of Biomolecular Engineering, Tagawa Laboratory, Japan).</i> <i>Die Daten wurden von Dr. Ahn (Tokyo Institute of Technology Graduate School of Bioscience and Biotechnology Department of Biomolecular Engineering, Tagawa Laboratory, Japan) zur Verfügung gestellt.</i></p>