

T4 Polynucleotide Kinase

Catalogue number: MB00801, 500 U

Description

T4 Polynucleotide Kinase catalyses the transfer of the terminal phosphate of ATP to 5'-hydroxyl termini of polynucleotides such as DNA and RNA, oligonucleotides and 3'-mononucleotides. In the presence of ADP, it will also catalyse the exchange of 5'-terminal phosphate groups and ATP. This enzyme also possesses a 3'-phosphatase and 2', 3' cyclic phosphodiesterase activity.

Storage

50 mM Tris-HCl, pH 7.6, 1.0 mM DTT, 0.1 mM EDTA, 50% (v/v) glycerol. Store at -20 °C.

Concentration: 30 U/μL

Unit definition

One enzyme unit incorporates 1 nmol of radiolabeled ATP into DNA substrate in 30 min at 37 °C under standard assay conditions.

Reaction buffer (10×): 0.5 M Tris-HCl, pH 7.6, 100 mM MgCl₂, 100 mM 2-mercaptoethanol. Supplied with dilution buffer: 50 mM Tris-HCl, pH 8.0.

Related products:

Product name	Cat. No.
Alkaline Phosphatase	MB018
Klenow Fragment of DNA polymerase I	MB009
Speedy Ligase	MB130
NZYGelpure	MB011

5' – End labelling protocol

1. This protocol is suitable for 10 pmol of oligonucleotide. Combine the following:

Oligonucleotide (10 pmol)	__ μL
[γ- ³² P]ATP*	__ μL
10× Reaction buffer	5 μL
Water	__ μL
T4 Polynucleotide Kinase (diluted)	2-5 units
Final volume	50 μL

2. Mix and incubate at 37 °C for 30 min.
 3. Terminate the reaction by heating at 65 °C for 10 minutes. Proceed with the separation of 5'-end-labeled oligonucleotide from precursor ATP by thin-layer or column chromatography according to standard methods.

* The number of pmol of [γ-³²P]ATP should be at least 2 times the number of pmol of 5'-ends.

Notes: Dilute the T4 Polynucleotide Kinase using 50 mM Tris-HCl, pH 8.0 (Dilution buffer). Adding concentrated enzyme can result in lower levels of incorporation. Do not store the enzyme in the diluted form. The labelling efficiency is influenced by the terminal 5' nucleotide. An oligonucleotide with G at the 5'-end labels about 6 fold higher than a 5'C and about 1.5 fold higher than a 5'T or 5'A.

Revised 12/12

Certificate of Analysis

Test	Result
Enzyme purity	Pass
Nucleases assays	Pass
Functional assay	Pass

Approved by:



José Prates
 Senior Manager, Quality Systems