

## GeneJET™ Gel Extraction Kit, #K0691, #K0692


**Note.** All steps should be carried out at room temperature. All centrifugations should be carried out in a microcentrifuge at  $\geq 12\,000 \times g$  (10 000-14 000 rpm, depending on the rotor type).

1



### Prepare the DNA for binding

Excise the gel slice with the DNA fragment and weigh.  
Add a 1:1 volume of **Binding Buffer** to the gel slice (volume/weight) (e.g. add 100  $\mu\text{l}$  of Binding buffer for every 100 mg of agarose gel).  
Incubate the gel mixture at 50-60°C for 10 min. Invert the tube every few min.



- Note:**
- For gels with an agarose content  $\geq 2\%$ , add 2:1 volumes of Binding Buffer to the gel slice.
  - If the DNA fragment is  $\leq 500$  bp or  $> 10$  kb, see complete protocol in the manual for additional recommendations.



2



### Bind DNA

Transfer the solution to the GeneJET™ Purification Column.  
Centrifuge for 30-60 s. Discard the flow-through.

- Note:** If the purified DNA will be used for sequencing see complete protocol in the manual for additional recommendations.



3



### Wash the column

Add 700  $\mu\text{l}$  of **Wash Buffer** and centrifuge for 1 min.  
Discard the flow-through.  
Centrifuge empty column for 1 min.



4



### Elute purified DNA

Place the column into a fresh 1.5 ml microfuge tube.  
Add 50  $\mu\text{l}$  of **Elution Buffer** to the column.  
Centrifuge for 1 min. Collect the flow-through.