

PRODUCT INFORMATION
Thermo Scientific
CloneJET PCR Cloning
Kit

#K1231 **20 rxns**
Lot **Expiry Date**
Store at -20 °C

www.thermoscientific.com/onebio

COMPONENTS OF THE KIT

| CloneJET PCR Cloning Kit | #K1231 20 cloning rxns | #K1232 40 cloning rxns |
|---|------------------------------|------------------------------|
| pJET1.2/blunt Cloning Vector (50 ng/μL) | 24 μL | 46 μL |
| 2X Reaction Buffer | 300 μL | 600 μL |
| T4 DNA Ligase (5 U/μL) | 24 μL | 46 μL |
| DNA Blunting Enzyme | 24 μL | 46 μL |
| pJET1.2 Forward Sequencing Primer, 10 μM aqueous solution | 50 μL | 100 μL |
| pJET1.2 Reverse Sequencing Primer, 10 μM aqueous solution | 50 μL | 100 μL |
| Control PCR Product (24 ng/μL) 976 bp, with 3'-dA overhangs | 8 μL | 12 μL |
| Water, nuclease-free | 1.25 mL | 1.25 mL |

STORAGE

All components of the Thermo Scientific™ CloneJET™ PCR Cloning Kit should be stored at -20°C.

CERTIFICATE OF ANALYSIS

All components of the kit were function tested in control experiment as described in the manual. A 2 μL aliquot of the ligation mixture was used to transform 50 μL of chemically competent XL1-Blue cells of >10⁸ cfu/μg DNA transformation efficiency. Cloning efficiency of the Control PCR Product into the pJET1.2/blunt was >2x10⁴ cfu/μg DNA. >90% of the recombinant plasmids contained the appropriate insert.

 Quality authorized by: Jurgita Zilinskiene

DESCRIPTION

The CloneJET PCR Cloning Kit is an advanced positive selection system for the highest efficiency cloning of PCR products generated with *Pfu* DNA polymerase, *Taq* DNA polymerase, Thermo Scientific™ DreamTaq™ DNA polymerase or other thermostable DNA polymerases. Additionally, any other DNA fragment, either blunt or sticky-end, can be successfully cloned using the kit. The kit features the novel positive selection cloning vector pJET1.2/blunt. This vector contains a lethal gene which is disrupted by ligation of a DNA insert into the cloning site. As a result, only cells with recombinant plasmids are able to propagate, eliminating the need for expensive blue/white screening. The vector contains an expanded multiple cloning site, as well as a T7 promoter for *in vitro* transcription. Sequencing primers are included for convenient sequencing of the cloned insert.

CLONING PRINCIPLE

pJET1.2/blunt is a linearized cloning vector, which accepts inserts from 6 bp to 10 kb. The 5'-ends of the vector contain phosphoryl groups, therefore, phosphorylation of the PCR primers is not required. Blunt-end PCR products generated by proofreading DNA polymerases can be directly ligated with the vector in just 5 min. PCR products with 3'-dA overhangs generated using *Taq* DNA polymerase or other non-proofreading thermostable DNA polymerase are blunted in 5 min with a proprietary thermostable DNA blunting enzyme (included in the kit) prior to ligation. All common laboratory *E. coli* strains can be directly transformed with the ligation product.

IMPORTANT NOTES

- Thoroughly mix every vial before use.
- The CloneJET PCR Cloning Kit is compatible with all PCR buffers supplied by Thermo Scientific.
- Gel-analyze the PCR product for specificity and yield before cloning.
- Specific PCR products of <1 kb appearing as one discrete band on the gel can be used for ligation directly from PCR reaction mixture without any purification.
- Do not use more than 1 μL of unpurified PCR product in the blunting or ligation reaction. Excess polymerase (*Taq*, *Pfu* or other) or salts from the PCR reaction mixture may result in background colonies and may reduce the efficiency of the cloning procedure.
- Gel purification of the PCR product (e.g. with Thermo Scientific™ GeneJET™ Gel Extraction Kit, #K0691) is recommended to increase the number of recombinants containing full length inserts in following cases:
 - PCR product is longer than 1 kb;
 - PCR product is contaminated with non-specific PCR products;
 - PCR product is contaminated with primer-dimers;
 - PCR template contains β-lactamase (ampicillin resistance) gene, which may result in background colonies on LB-ampicillin agar plates.
- For efficient cloning of gel-purified DNA fragments, it is important to avoid DNA damage by ethidium bromide and UV light. Use a long wavelength UV (360 nm) light-box when excising DNA from the agarose gel. When using a short-wavelength (254-312 nm) light-box, limit DNA exposure to UV to a few seconds. Keep the gel on a glass plate or on a plastic plate during UV illumination. Alternatively, use dyes, like crystal violet, to visualize DNA in ambient light (1, 2).
- The kit performs well over a wide range of insert/vector molar ratios (0.5:1 to 15:1). However the optimal insert/vector ratio is 3:1. Vector pJET1.2/blunt is supplied at a concentration of 0.05 pmol DNA ends/μL. To calculate optimal amount of the PCR product for ligation (0.15 pmol of DNA ends respectively), refer to Table 1 or use dedicated software (e.g., www.thermoscientific.com/reviewer) for calculations.
- For PCR products >3 kb, ligation can be prolonged to 30 min. Ligation times longer than 30 min are not recommended and may decrease cloning efficiency.
- Success in cloning of long PCR products may also depend on the DNA sequence of the insert. PCR products may contain toxic

sequences not tolerated by *E.coli*, therefore multicopy vectors like pJET1.2 may not be suitable cloning these PCR products.

Table 1. Recommended amount of PCR product for the ligation reaction

| Length of PCR product (bp) | Optimal PCR product quantity for ligation reaction, (0.15 pmol ends) |
|----------------------------|--|
| 100 | 5 ng |
| 300 | 15 ng |
| 500 | 25 ng |
| 1000 | 50 ng |
| 2000 | 100 ng |
| 3000 | 150 ng |
| 4000 | 200 ng |
| 5000 | 250 ng |

CLONING PROTOCOLS

Blunt-End Cloning Protocol

- For cloning blunt-end PCR products generated by proofreading DNA polymerases, such as *Pfu* DNA polymerase. If the DNA end structure of the PCR products is not specified by the supplier of the DNA polymerase, follow the Sticky-End Cloning Protocol.
- For cloning of blunt-end DNA fragments generated by restriction enzyme digestion. Gel-purify the DNA fragment prior to ligation and use in a 3:1 molar ratio with pJET1.2/blunt (see Table 1).

1. Set up the ligation reaction **on ice**:

| Component | Volume |
|---|------------------------|
| 2X Reaction Buffer | 10 μL |
| Non-purified PCR product or purified PCR product/other blunt-end DNA fragment | 1 μL 0.15 pmol ends |
| pJET1.2/blunt Cloning Vector (50 ng/μL) | 1 μL (0.05 pmol ends) |
| Water, nuclease-free | up to 19 μL |
| T4 DNA Ligase | 1 μL |
| Total volume | 20 μL |

Vortex briefly and centrifuge for 3-5 s.

2. Incubate the ligation mixture at room temperature (22°C) for 5 min. **Note.** For PCR products >3 kb, ligation can be prolonged to 30 min.
3. Use the ligation mixture directly for transformation **Note.** Keep the ligation mixture at -20°C if transformation is postponed. Thaw on ice and mix carefully before transformation.

Sticky-End Cloning Protocol

- For cloning PCR products with 3'-dA overhangs generated by *Taq* DNA polymerase, DreamTaq DNA polymerase or enzyme mixtures containing *Taq* DNA polymerase.
- For cloning PCR products when DNA end structure of the generated PCR products is not specified by the supplier of the DNA polymerase.
- For cloning DNA fragments with 5'- or 3'-overhangs generated by restriction enzyme digestion. Gel-purify the DNA fragment prior to ligation and use in a 3:1 molar ratio with pJET1.2/blunt (see Table1). **Note.** The DNA Blunting Enzyme is a proprietary thermostable DNA polymerase with proofreading activity. It will remove 3'-overhangs and fill-in 5'-overhangs. Nucleotides for the blunting reaction are included in the reaction buffer.

1. Set up the blunting reaction **on ice**:

| Component | Volume |
|--|--------------------------|
| 2X Reaction Buffer | 10 μL |
| Non-purified PCR product or purified PCR product/other sticky-end DNA fragment | 1 μL (0.15 pmol ends) |
| Water, nuclease-free | to 17 μL |
| DNA Blunting Enzyme | 1 μL |
| Total volume | 18 μL |

Vortex briefly and centrifuge for 3-5 s.

2. Incubate the mixture at 70 °C for 5 min. Chill on ice.
3. Set up the ligation reaction **on ice**. Add the following to the blunting reaction mixture:

| Component | Volume |
|---|--------------------------|
| pJET1.2/blunt Cloning Vector (50 ng/μL) | 1 μL (0.05 pmol ends) |
| T4 DNA Ligase | 1 μL |
| Total volume | 20 μL |

Vortex briefly and centrifuge for 3-5 s to collect drops.

4. Incubate the ligation mixture at room temperature (22 °C) for 5 min. **Note.** For PCR products >3 kb, ligation can be prolonged to 30 min.
5. Use the ligation mixture directly for transformation **Note.** Keep the ligation mixture at -20 °C if transformation is postponed. Thaw on ice and mix carefully before transformation.

Transformation

- The CloneJET PCR Cloning Kit is compatible with all common *E.coli* laboratory strains. Transformation of competent *E.coli* cells with the ligation mixture can be performed using different transformation methods (Table 2 and 3).
- The number of transformants on the plates directly depends on the transformation efficiency of the competent cells.
- For successful cloning, competent *E.coli* cells should have an efficiency of at least 1 × 10⁶ cfu/μg supercoiled plasmid DNA. To check the efficiency, prepare a control transformation with 0.1 ng of a supercoiled vector DNA, e.g., pUC19 DNA, #SD0061 (Table 2).
- For transformation of the ligation mixture, refer to Table 3.

Table 2. Evaluation of transformation efficiency of competent cells (control transformation)

| Transformation method | Number of transformants per μg of supercoiled plasmid DNA | Amount of pUC19 DNA to yield ~ 1000 colonies per plate | Volume of competent cells |
|--|---|--|---------------------------|
| Thermo Scientific™ TransformAid™ Bacterial Transformation Kit, #K2710* | ~1 × 10 ⁷ | 0.1 ng | 50 μL |
| Calcium chloride transformation | ~1 × 10 ⁶ | 1 ng | 50 μL |
| Electro-transformation | ~1 × 10 ⁹ | 0.01 ng | 40 μL |

* XL1-Blue, ER2267, ER1727 *E. coli* strains are the best strains for transformation with TransformAid™ Bacterial Transformation Kit. DH10B, DH5α and TOP10 strains are not efficient with TransformAid but are recommended for calcium chloride transformation or electroporation.

Table 3. Recommendations for transformation of ligation mixture

| Transformation method | Treatment of the ligation mixture before transformation | Volume of the ligation mixture for transformation | Volume of competent cells for transformation |
|---|--|---|--|
| TransformAid Bacterial Transformation Kit, #K2710 | Do not heat ligation mixture prior used for transformation. | ≤2.5 µL | 50 µL |
| Calcium Chloride Transformation | Do not heat ligation mixture prior used for transformation. | ≤5 µL | 50 µL |
| Electro-transformation | Spin column or chloroform extraction (see protocol below) Do not heat ligation mixture prior used for transformation. | 0.5-1 µL of purified ligation mixture | 40 µL |

Analysis of recombinant clones

Analyze 4-6 colonies for the presence and orientation of the DNA insert using one of the following methods:

Colony PCR

Use the following protocol for colony screening by PCR if the cloned PCR fragment is shorter than 3 kb. For longer inserts, perform restriction analysis.

1. Prepare enough PCR master mix for the number of colonies analyzed plus one extra. For each 20 µL reaction, mix the following reagents:

| Component | Using Taq DNA Polymerase | Using 2X PCR Master Mix |
|--|--------------------------|-------------------------|
| 10X Taq buffer | 2.0 µL | – |
| dNTP mix, 2 mM each | 2.0 µL | – |
| 25 mM MgCl ₂ | 1.2 µL | – |
| pJET1.2 Forward Sequencing Primer, 10 µM | 0.4 µL | 0.4 µL |
| pJET1.2 Reverse Sequencing Primer, 10 µM | 0.4 µL | 0.4 µL |
| Water, nuclease-free | 13.9 µL | 9.2 µL |
| Taq DNA Polymerase 5 u/µL, #EP0401 or DreamTaq Green DNA Polymerase, #EP0711 | 0.1 µL | – |
| PCR Master Mix (2X), #K0171 or DreamTaq Green PCR Master Mix (2X), #K1081 | – | 10 µL |
| Total volume | 20 µL | 20 µL |

2. Mix well. Aliquot 20 µL of the mix into the PCR tubes on ice.
3. Pick a small amount of individual colony and resuspend in 20 µL of the PCR master mix. *Picking a colony, take care not to touch the agar.*

Note. Due to considerable amount of recircularised vector plated on the surface of plate, colony PCR may give some false-negative results. Prior to clone analysis propagate short strikes of individual colonies on ampicillin plates. Then use small amount of each for colony PCR.

4. Perform PCR: 95 °C, 3 min; 94 °C, 30 s, 60 °C, 30 s, 72 °C 1 min/kb; 25 cycles.
5. Analyze on an agarose gel for the presence of the PCR product.

Restriction analysis

Isolate plasmid DNA from an overnight bacterial culture. To speed up the process and to assure the quality of purified plasmid DNA, use the Thermo Scientific™ GeneJET™ Plasmid Miniprep Kit, #K0503. To digest DNA from recombinant clones in just 5 min, use Thermo Scientific™ FastDigest™ restriction enzymes.

Sequencing

Use the pJET1.2 Forward Sequencing Primer or pJET1.2 Reverse Sequencing Primer supplied with the kit to sequence the cloned insert.

CONTROL CLONING EXPERIMENT

The control reaction should be used to verify the efficiency of the blunting and ligation steps. The 976 bp control PCR product (nucleotide sequence is available at www.thermoscientific.com/onebio) has been generated with Taq DNA polymerase, which adds extra nucleotides to the 3'-end. Therefore, the **Sticky-End Protocol** must be followed.

1. Set up the blunting reaction on ice:

| Component | Volume |
|--------------------------------|--------|
| 2X Reaction Buffer | 10 µL |
| Control PCR Product (24 ng/µL) | 2 µL |
| Water, nuclease-free | 5 µL |
| DNA Blunting Enzyme | 1 µL |
| Total volume | 18 µL |

- Vortex briefly and centrifuge for 3-5 s to collect drops.
2. Incubate the mixture at 70 °C for 5 min. Chill on ice.
3. Set up the ligation reaction on ice. Add the following to the blunting reaction mixture:

| Component | Volume |
|---|--------|
| pJET1.2/blunt Cloning Vector (50 ng/µL) | 1 µL |
| T4 DNA Ligase | 1 µL |
| Total volume | 20 µL |

- Vortex briefly and centrifuge for 3-5 s to collect drops.
4. Incubate the ligation mixture at room temperature (22 °C) for 5 min.
5. Use the ligation mixture directly for transformation. Keep the ligation mixture at -20 °C if transformation is postponed. Thaw on ice and mix carefully before transformation.

Analyze colonies by colony PCR. At least 9 of 10 analyzed colonies should contain recombinant plasmid with the 976 bp insert. The number of transformants depends on the transformation efficiency of the *E. coli* cells. Verify the transformation efficiency by transforming supercoiled plasmid, e.g., pUC19 DNA, #SD0061 in parallel. Refer to Table 2 for correct control transformations.

MAP AND FEATURES OF pJET1.2/blunt CLONING VECTOR

The pJET1.2/blunt cloning vector has been linearized with Eco321 (EcoRV) (GenBank/EMBL Accession number EF694056). The blunt ends of the vector contain 5'-phosphoryl groups. The nucleotide sequence of pJET1.2/blunt is available at www.thermoscientific.com/onebio. For complete list of enzymes that cut or do not cut the pJET1.2/blunt see www.thermoscientific.com/reviewer.




Fig. 1. pJET1.2/blunt Vector Map.
Genetic elements of the pJET1.2/blunt cloning vector

| Element | Function | Position (bp) |
|-----------------------------------|--|---------------|
| rep (pMB1) | Replicon (rep) from the pMB1 plasmid responsible for the replication of pJET1.2 | 1762-1148 |
| Replication start | Initiation of replication | 1162±1 |
| bla (Amp ^R) | β-lactamase gene conferring resistance to ampicillin. Used for selection and maintenance of recombinant <i>E. coli</i> cells | 2782-1922 |
| eco47IR | Lethal gene <i>eco47IR</i> enables positive selection of recombinant plasmid | 753-16 |
| P _{lacUV5} | Modified P _{lac} promoter for expression of the <i>eco47IR</i> gene at a level sufficient to allow for positive selection | 892-769 |
| T7 promoter | T7 RNA polymerase promoter for <i>in vitro</i> transcription of the cloned insert | 305-324 |
| Multiple cloning site (MCS) | Mapping, screening and excision of the cloned insert | 422-328 |
| Insertion site | Blunt DNA ends for ligation with insert | 371-372 |
| Primer binding sites: | | |
| pJET1.2 forward sequencing primer | Sequencing of insert, colony PCR. Sequence: 5'-CGACTCACTATAGGGAGAGCGGC-3' | 310-332 |
| pJET1.2 reverse sequencing primer | Sequencing of insert, colony PCR. Sequence: 5'-AAGAACATCGATTTTCATGGCAG-3' | 428-405 |

RECIPES & SUPPLEMENTARY PROTOCOLS**Chloroform extraction of the ligation mixture prior electroporation**

1. Add an equal volume (20 µL) of the chloroform to the ligation mixture. Mix well.
2. Centrifuge at 10,000 rpm for 3 min at room temperature.
3. Carefully transfer the upper aqueous phase to a fresh tube.

4. Use 1 µL of the purified mixture for transformation of 40 µL of electrocompetent cells.
- Keep the purified ligation mixture at -20 °C if transformation is postponed. Thaw on ice and mix carefully before transformation.

Ampicillin stock solution (50 mg/mL)

1. Dissolve 2.5 g ampicillin sodium salt in 50 mL of deionized water.
2. Filter sterilize and store in aliquots at 4 °C.

LB-ampicillin plates

1. Prepare LB-agar Medium (1 liter), weigh out:
Bacto® Tryptone 10 g,
Bacto Yeast extract 5 g,
NaCl 5 g.
2. Dissolve in 800 mL of water, adjust pH to 7.0 with NaOH and add water to 1000 mL.
3. Add 15 g of agar and autoclave.
4. Allow the medium to cool to 55 °C.
5. Add 2 mL of ampicillin stock solution (50 mg/mL) to a final concentration of 100 µg/mL.
6. Mix gently and pour plates.

TROUBLESHOOTING

For troubleshooting see the complete version of product information for CloneJET PCR Cloning Kit at www.thermoscientific.com/onebio.

References

1. Rand, K.N., Crystal Violet can be used to Visualize DNA Bands during Gel Electrophoresis and to Improve Cloning Efficiency, Elsevier Trends Journals Technical Tips, Online, T40022, 1996.
2. Adkins, S., Burmeister, M., Visualization of DNA in agarose gels and educational demonstrations, Anal Biochem., 240 (1), 17-23, 1996.

NOTICE TO PURCHASER:

Use of this product is covered by US Patent No. 6,127,155. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim, no right to perform any patented method and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting outlicensing@lifetech.com or Out Licensing, Life Technologies Inc., 5791 Van Allen Way, Carlsbad, California 92008.

LIMITED USE LABEL LICENSE: Internal Research and Development Use Only.

The purchase of this product conveys to the buyer the limited, non-exclusive, non-transferable right (without the right to resell, repackage, or further sublicense) to use this product for internal research and development purposes. No other license is granted to the buyer whether expressly, by implication, by estoppel or otherwise. In particular, the purchase of the product does not include or carry any right or license to use, develop, or otherwise exploit this product commercially and no rights are conveyed to the buyer to use the product or components of the product for purposes including but not limited to provision of services to a third party, generation of commercial databases or clinical diagnostics. This product is sold pursuant to authorization from Thermo Fisher Scientific and Thermo Fisher Scientific reserves all other rights. For information on purchasing a license for uses other than internal research and development purposes, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies Inc., 5791 Van Allen Way, Carlsbad, California 92008.

PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively for research purposes and *in vitro* use only. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals. Please refer to www.thermoscientific.com/onebio for Material Safety Data Sheet of the product.