

## NZYTaq II DNA polymerase

**Catalogue number:** MB35401, 500 U  
 MB35402, 1000 U  
 MB35403, 2500 U

### Description

NZYTaq II DNA polymerase belongs to a new generation of *Taq*-derived DNA polymerases optimized for standard PCR applications. The enzyme was engineered to produce high DNA yields in shorter PCR running times (15-30 s/kb extension) under minimal optimization conditions. NZYTaq II DNA polymerase lacks 3'→5' exonuclease activity and supports the robust and reliable amplification of a wide range of DNA templates up to 6 kb. The enzyme was optimized to provide higher sensitivity, allowing amplification of different DNA fragments from as little as 5 pg of human genomic DNA. Resulting PCR products have an A-overhang and are suitable for cloning with NZYTaq's TA PCR cloning kits (MB053 or MB137).

### Storage temperature

NZYTaq II DNA polymerase should be stored at -20 °C in a constant temperature freezer. NZYTaq II DNA polymerase will remain stable till the expiry date if stored as specified.

### Unit definition

One unit is defined as the amount of enzyme required to catalyse the incorporation of 10 nmoles of dNTPs into acid insoluble material in 30 minutes at 72 °C.

**Enzyme concentration:** 5 U/μL

**Magnesium Chloride solution:** the provided 50 mM MgCl<sub>2</sub> solution allows users to optimize Mg<sup>2+</sup> concentration in different PCR set ups. In general, NZYTaq II DNA polymerase works effectively with a 2.5 mM MgCl<sub>2</sub> concentration. Vortex the MgCl<sub>2</sub> solution thoroughly after thawing.

### Standard Protocol

The following standard protocol serves as a general guideline and a starting point for any PCR amplification. Optimal reaction conditions (e.g. concentration of DNA polymerase, primers, MgCl<sub>2</sub> and template DNA) may vary, although PCR optimization is usually not required. In case you need to fine-tune PCR conditions, recommended variations of each PCR component are provided in brackets in the table below.

1. Gently mix and briefly centrifuge all components after thawing. On ice, in a sterile, nuclease-free microcentrifuge

tube, prepare a mixture for the appropriate number of PCR reactions. Add water first and the remaining components in the order specified in the table below. A single reaction mixture of 50 μL should combine the following components:

Reaction buffer, 10× (provided)	5 μL
MgCl <sub>2</sub> , 50 mM (provided)	2.5 (1.5-4.0) mM
dNTPs mix	0.4 (0.25-0.5) mM
Primers (see below)	0.25 (0.1-0.5) μM
Template DNA (see below)	5 pg-0.5 μg
NZYTaq II (5 U/μL)	0.25-1 μL
Nuclease-free water	up to 50 μL

2. Mix and quickly pulse the reactions.

3. Perform PCR using the following cycling parameters:

Cycle step	Temp.	Time	Cycles
Initial denaturation	95 °C	3 min	1
Denaturation	94 °C	30 s	25-35
Annealing	*	30 s	
Extension	72 °C	15-30 s/kb	
Final Extension	72 °C	5-10 min	1

\*Annealing temperature should be optimized for each primer set based on the primer T<sub>m</sub>; typically it should be T<sub>m</sub>-5 °C.

4. Analyse PCR products through agarose gel electrophoresis (0.7-1.2%, w/v) and visualise with GreenSafe Premium (MB132) or any other means.

### Primer Design

PCR primers generally range in length from 15–30 bases and are designed to flank the region of interest. Primers should contain 40–60% GC, and care should be taken to avoid sequences that might produce internal secondary structure. The 3'-ends of the primers should not be complementary to avoid the production of primer-dimers. Primer-dimers unnecessarily remove primers from the reaction and result in an unwanted polymerase reaction that competes with the desired reaction. Avoid three G or C nucleotides in a row near the 3'-end of the primer, as this may result in non-specific primer annealing. Ideally, both primers should have nearly identical melting temperatures (T<sub>m</sub>), allowing their annealing with the denatured template DNA at roughly the same temperature.

### DNA template

The optimal amount of starting material may vary depending on the quality and complexity of template DNA. In general, we recommend using 20-500 ng of genomic DNA templates, although the enzyme is sensitive enough to amplify DNA fragments from as little as 5 pg of human gDNA, for example. Lower amounts of DNA template (typically 1-20 ng) may be used for amplification of lambda or plasmid DNA or even 5-20 ng for amplification of multicopy chromosomal

genes. When using a cDNA synthesis reaction as template do not exceed 10% of the final PCR reaction volume.

## Quality control assays

### Purity

NZYTaq II DNA polymerase purity is >90% as judged by SDS polyacrylamide gel electrophoresis followed by Coomassie Blue staining.

### Genomic DNA contamination

NZYTaq II DNA polymerase must be free of any detectable genomic DNA contamination as evaluated through PCR.

### Nuclease assays

0.2-0.3 µg of pNZY28 plasmid DNA are incubated with 5 U of NZYTaq II DNA polymerase, in 1× Reaction Buffer, for 14-16 hours at 37 °C. Following incubation, the DNA is visualised on a GreenSafe Premium-stained agarose gel. There must be no visible nicking or cutting of the nucleic acid. Similar tests are performed with NZYTaq II reaction buffer and MgCl<sub>2</sub> solution.

### Functional assay

NZYTaq II DNA polymerase is extensively tested for performance in a polymerase chain reaction (PCR) of different-sized DNA fragments (1 and 2.5 kb) from human genomic DNA in the presence of 10× Reaction Buffer and MgCl<sub>2</sub> solution. The resulting PCR products are visualized as single bands in a GreenSafe Premium-stained agarose gel.

## Troubleshooting

### No product amplification or low yield

- Inadequate annealing temperature

The reaction mix composition may affect the melting properties of primers and DNA. Adjust the annealing temperature to accommodate the primer with the lowest melting temperature (5 °C to 10 °C lower than T<sub>m</sub>).

- Presence of PCR inhibitors

Some DNA isolation procedures, particularly genomic DNA isolation, can result in the co-purification of PCR inhibitors. Reduce the volume of template DNA in reaction or dilute template DNA prior to adding to the reaction. Diluting samples even 1:10,000 has been shown to be effective in improving results, depending on initial DNA concentration.

- Additives required

Adding PCR-enhancing agents (NZYTaq 5× Optimizer Solution – MB060 or NZYTaq 2× GC-Enhancer Solution – MB143) may improve yield while allowing the amplification of difficult templates.

### Presence of non-specific bands

- Non-specific annealing of primers

Adjust annealing conditions and/or design another set of primers, by increasing the length and avoiding complementary sequences.

- Mg<sup>2+</sup> concentration is too high

Generally, 2-3 mM MgCl<sub>2</sub>, typically 2.5 mM final concentration, works well for the majority of PCR reactions. Optimal concentration depends on target template, buffer and dNTPs. Optimize magnesium concentration by supplementing MgCl<sub>2</sub> in 0.5 increments up to 4 mM.

Revised 06/17

## Certificate of Analysis

Test	Result
Enzyme purity*	Pass
Genomic DNA contamination*	Pass
DNase contamination	Pass
Functional assay	Pass

\*These assays were performed exclusively with the enzyme

Approved by:



José Prates  
Senior Manager, Quality Systems